




Improving Department Results with Single Exposure Dual-Energy Subtraction X-Ray: Clinical, Operational and Financial Benefits


KA Imaging
July 27, 2021

Karim S Karim PhD PEng MBA FEIC
Founder and Chief Technology Officer

Jay Potipcoe, MRT(R), BA
Applications Specialist



Dual-Energy (DE) X-Ray Subtraction (1970s)



Studies have shown that DE enables accurate diagnosis by removing bones/soft tissue from X-ray image.



LOW ENERGY



HIGH ENERGY



SOFT



HARD



Problems with Dual Energy

Access

- Higher radiation per scan → Limited use cases
- **Not portable**
- Cannot substitute standard chest X-ray


Quality

- DE images unreadable up to 20% of the time due to misregistration errors (**motion artifact**)

Cost

- Fixed system available from very few manufacturers
- **Large upfront costs** (acquisition, installation and maintenance)

DE solutions to date have seen slow adoption because of motion artifacts, higher radiation, change in workflow and expensive X-ray systems



Dual Energy Motion Artifact
RSNA Radiographics(2006) vol. 26 pp. 79-91



Dual-Energy Solutions: Dual-layer Detector Technology (1990's)



X-RAY
SOURCE



Conventional
X-Ray source

Energy Separation performed in the detector

X-RAY
DETECTOR



✖

**DOSE
EFFICIENCY**
 $DQE(u, v)$

POOR

- Significant loss of photons in mid-filter

✖

**DUAL-
ENERGY
IMAGE
QUALITY**

POOR

- Tradeoff between mid-filter thickness, dose efficiency and spectral separation



Jay Potipcoe, BA, MRT



Dr. Karim Karim, Ph.D., PEng, MBA

Dual-Energy Solutions: Bone Suppression Software (2010's)



X-RAY
SOURCE



Conventional
X-Ray source

Morphology based software solution (not true dual-energy) used to perform bone and tissue separation

X-RAY
DETECTOR



✓

GOOD

Same as chest X-ray

✗

POOR

- Limited use cases (i.e. nodules)
- Limited views (i.e. PA/AP only)
- Cannot effectively distinguish between soft tissue and bone



Jay Potipcoe, BA, MRT



Dr. Karim Karim, Ph.D., PEng, MBA

Dual-Energy Solutions: Triple-layer Detector Technology (2020's)

X-RAY
SOURCE




Conventional
X-Ray source

X-RAY
DETECTOR

Energy Separation performed in
the detector




Single X-ray Exposure



- EXCELLENT**
- Three sensor layers capture more photons
 - Better DQE than most chest X-rays
- EXCELLENT**
- Layers optimized for high spectral separation
 - Zero motion artifact



How Do Three Layers Work?



- Single 120 kVp exposure generates three distinct spectral X-ray images
- Thickness for each layer can be optimized for best spectral separation without losing X-ray dose efficiency
- X-ray signals from all three layers are added to get the best-in-class dose efficiency X-ray image
- Signals from Layer 1 and Layer 3 are subtracted to yield dual energy images
- Simultaneous acquisition of three images results in zero motion artifact




Single Exposure Dual Energy Solution



Three images from one chest X-ray exposure maintains same clinical technique and dose as digital chest X-ray



KA
Standard X-Ray



High Energy Bone X-Ray



Low Energy Soft Tissue X-Ray



Comparison Between Different Solutions



	Single Exposure Dual-Energy	X-ray + Bone Suppression Software	Dual Exposure Dual Energy	Digital Chest X-ray
Clinical Use Cases	More than Chest X-ray	Limited to nodules	Less than Chest X-ray	Standard
Number of X-rays (Radiation)	Single	Single	Double	Single
Use of Grid	As per clinical technique	As per clinical technique	Always required	As per clinical technique
Better Sensitivity than Chest X-ray	Yes	Marginal	Yes	No
Distinguish Calcium from Soft Tissue	Yes (physics solution)	Marginal (software solution)	Yes (physics solution)	No
AP, PA, Lateral and Oblique Views	Yes	No	No	Yes
Portable and Mobile Applications	Yes	Yes	No	Yes
Motion Artifacts	Absent	Absent	Present	N/A





Motion Artifact Comparison



Single exposure DE produces DE images free of motion artifacts

**Solution 1:
Dual exposure DE**




**Solution 2:
Dual exposure DE**



**Solution 3:
Dual exposure DE**



**Solution 4:
Single exposure DE**








CASE A: Hidden Apical Lesion in RUL



Conventional DR X-Ray

Soft Tissue Image

Bone Image



-This DR image was read as normal by the radiologist, but the patient actually has an apical lung lesion that was missed, due to overlapping bone. When shown the DE images, pathology was seen by the radiologist. Without DE, this would have gone unnoticed, losing valuable proactive treatment time.

-Incidentally, an older rib fracture also became visible on right side.



Jay Potipcoe, BA, MRT



Dr. Karim Karim, Ph.D., PEng, MBA




CASE B: Hidden Mass in RLL Obscured by Degenerative Bone Disease (PA)



Conventional DR X-Ray

Soft Tissue Image

Bone Image



-This granuloma in the DR image can be classified as being calcified, due to its appearance in the Bone Image.

-The PortoCath line is very easy to see on the Bone Image

-There is, however, a mass behind the heart not seen in this PA image. Because KA Imaging's Single Exposure DE detector can perform lateral chest imaging, the mass will be seen in that view.



CASE B: Hidden Mass in RLL Obscured by Degenerative Bone Disease (Lat)



Conventional DR X-Ray

Soft Tissue Image

Bone Image



-The mass seen behind the heart is superimposed over a calcified joint in the conventional DR image. However, the Soft Tissue image clearly displays the mass to the radiologist, by removing the overlapping bone.

-The Portocath line and sternum are also nicely differentiated in the Bone Image.




CASE C: Focal Opacity in RUL



Conventional DR X-Ray



Soft Tissue Image



CT Image

-This case involved a 51-year-old female leukemia patient. The Conventional DR X-Ray was given to the radiologist and read as normal. However, upon subsequently being shown the DE Soft Tissue image, the radiologist noticed a highlighted focal opacity indicating pneumonia (confirmed on CT). This patient would not have received prompt treatment for the pneumonia, if not for the DE image.




CASE C: Focal Opacity in RUL



Conventional DR X-Ray




Soft Tissue Image



CT Image




-This case involved a 51-year-old female leukemia patient. The Conventional DR X-Ray was given to the radiologist and read as normal. However, upon subsequently being shown the DE Soft Tissue image, the radiologist noticed a highlighted focal opacity indicating pneumonia (confirmed on CT). This patient would not have received prompt treatment for the pneumonia, if not for the DE image.




CASE D: Left Perihilar Opacity



Conventional DR X-Ray



Soft Tissue Image



CT Image


-The Conventional DR image of this 74-year-old male with leukemia was read as showing a possible fungal pneumonia in the left hilar region. The radiologist placed confidence in diagnosis at a mere 1/5. Given the DE Soft Tissue image, the radiologist placed confidence as 5/5 for left hilar fungal pneumonia. Two other bilateral opacities were also noted in the Soft Tissue image (RUL and LLL), that were not shown in the Conventional DR image. Diagnosis was confirmed by CT.




CASE D: Left Perihilar Opacity




Conventional DR X-Ray



Soft Tissue Image



CT Image



-The Conventional DR image of this 74-year-old male with leukemia was read as showing a possible fungal pneumonia in the left hilar region. The radiologist placed confidence in diagnosis at a mere 1/5. Given the DE Soft Tissue image, the radiologist placed confidence as 5/5 for left hilar fungal pneumonia. Two other bilateral opacities were also noted in the Soft Tissue image (RUL and LLL), that were not shown in the Conventional DR image. Diagnosis was confirmed by CT.



Jay Potipcoe, BA, MRT



Dr. Karim Karim, Ph.D., PEng, MBA


CASE E: Coronary Stents



Conventional DR X-Ray



Bone Image



-Cardiac calcifications and stents can be seen clearly using a DE Bone Image in lateral configuration





CASE F: Coronary Calcium



Conventional DR X-Ray

Soft Tissue Image

Bone Image



- Calcified Coronary vessels clearly seen in lateral image
- Clearly seen PortoCath line, sternum, and sternal wire ties in Bone Image




CASE G: Pneumothorax




Conventional DR X-Ray



Soft Tissue Image



Bone Image



-This 26-year-old male presented to the ER with difficulty breathing and chest pain. A single exposure dual energy X-ray detector was used to generate all 3 images, and clearly displays the left apical pneumothorax, best seen in the Soft Tissue Image.




Benefits of Dual Energy Subtraction – Pneumonia



		Sensitivity		Specificity	
		Chest X-ray	Dual Energy	Chest X-ray	Dual Energy
Interstitial disease	Lung emphysema	45%	75%	84%	83%
	Scarring	77%	96%	76%	70%
	Reticular lung changes	81%	92%	87%	87%
Infective consolidation		41%	61%	98%	97%
Lung atelectasis		22%	67%	99%	99%

Martini, Katharina, Marco Baessler, Stephan Baumüller, and Thomas Frauenfelder. "Diagnostic accuracy and added value of dual-energy subtraction radiography compared to standard conventional radiography using computed tomography as standard of reference." PloS one 12, no. 3 (2017): e0174285.






Case Study #2: Indeterminate Nodules

Solitary Pulmonary Nodules (SPN) are prevalent in **0.15%** of the general population

150 million chest X-Rays performed in the US in 2018



Earlier and Easier Detection = Better Patient Outcomes, More Downstream Revenue, and Enhanced Risk Avoidance



Jay Potipcoe, BA, MRT



Dr. Karim Karim, Ph.D., PEng, MBA

Benefits of Dual Energy Subtraction - Nodules



	Sensitivity (%)		Specificity		Accuracy		PPV		NPV	
	Chest X-ray	Dual Energy	Chest X-ray	Dual Energy	Chest X-ray	Dual Energy	Chest X-ray	Dual Energy	Chest X-ray	Dual Energy
Small lung nodules < 20 mm, Solid, Part-solid, Non-solid	47.6	63	72.5	72.5	56.7	66.5	75.6	80.1	44.4	53.7

Oda, Seitaro, Kazuo Awai, Kenji Suzuki, Yumi Yanaga, Yoshinori Funama, Heber MacMahon, and Yasuyuki Yamashita. "Performance of radiologists in detection of small pulmonary nodules on chest radiographs: effect of rib suppression with a massive-training artificial neural network." American Journal of Roentgenology 193, no. 5 (2009): W397-W402.

Dual energy increases sensitivity by 33% relative to X-ray



Early detection for proactively finding cases



Replacement for Chest X-ray in Outpatient Settings



Jay Potipcoe, BA, MRT




Dr. Karim Karim, Ph.D., PEng, MBA

Click to save a picture to your desktop.

Case Study #3: Small-Volume Pneumothorax

Affects up to **35** out of every 100,000 people

Cost of treatment is ~ **\$25,000**/case.
Cost goes up if case worsens, due to misdiagnosis





Benefits of Dual Energy Subtraction - Pneumothorax



	Chest X-ray	Dual Energy
	# of cases found in study	# of cases found in study
Interns	6.2	8.8
Residents	9.1	11
Senior Physicians	9.5	10.9
Non-radiologists	9.2	10.6
Radiologists	9.5	10.6

Urbaneja, Ayla, Gauthier Dodin, Gabriela Hossu, Omar Bakour, Rachid Kechidi, Pedro Gondim Teixeira, and Alain Blum. "Added value of bone subtraction in dual-energy digital radiography in the detection of pneumothorax: impact of reader expertise and medical specialty." Academic radiology 25, no. 1 (2018): 82-87.

Dual energy detected up to 30% more small volume pneumothoraces



Early detection leads to cost savings and better patient outcomes



Replaces Portable X-ray in Critical Care Units



Clinical Use Cases with High Economic Impact



Use Case	Issue	Extra Cases Found Relative to Chest X-ray over 5 years	Improved by DE
Solitary Pulmonary Nodule	Lung is #1 Killer cancer worldwide	60 (US data)	Increase in sensitivity by 33% for nodules < 20mm
Coronary Calcium	Cardiac disease is the #1 Killer in developed countries	205 (US data)	Presence of calcium in DE bone image is strongly correlated with calcium in CT
Pneumonia	33% of patients misdiagnosed, CXR not sensitive	90 (US data)	Sensitivity increases for infectious consolidation by 50% and atelectasis by 3X
Pneumothorax	>50,000 due to errors in the USA, 20% die	15 (US data)	Up to 30% more small volume pneumothoraces found even by residents and interns
Tuberculosis	1.5 million die annually in LMICs	> 300 (in a high TB prevalence setting in Asia)	Up to 30% increase in sensitivity and NPV for consolidation and cavitation






Jay Potipcoe, BA, MRT



Dr. Karim Karim, Ph.D., PEng, MBA

Click to save a picture to your desktop.

AI vs. Dual Energy



CASE REPORT			
Abnormality Score: 61%			
Atel	Atelectasis	48%	right upper
Calc	Calcification	Low	
Con	Cardiomegaly	56%	present
Con	Consolidation	21%	left lower
Fib	Fibrosis	61%	right upper
MIW	Mediastinal Widening	Low	
Nod	Nodule	56%	right upper
PEF	Pleural effusion	Low	
Ppne	Pneumoperitoneum	Low	
Pne	Pneumothorax	Low	

CASE REPORT			
Abnormality Score: 61%			
Atel	Atelectasis	48%	right upper
Calc	Calcification	Low	
Con	Cardiomegaly	56%	present
Con	Consolidation	21%	left lower
Fib	Fibrosis	61%	right upper
MIW	Mediastinal Widening	Low	
Nod	Nodule	56%	right upper
PEF	Pleural effusion	Low	
Ppne	Pneumoperitoneum	Low	
Pne	Pneumothorax	Low	



This is an example of a commercial AI chest product.

The image on the left is a DR image and the two images on the right are the outputs of the AI algorithms.



Compared with Single Exposure Dual Energy



AI missed nodules in the lower part of the image that are near the mediastinum

There is also pneumonia seen in the soft tissue image that was not identified by the AI



Looking to the Future: Coronary Calcium Scoring (New Use Case for X-ray)




-Data shows dual energy X-ray is adept at displaying coronary calcifications. There is exciting potential for using DE X-ray for calcium scoring^[1]

-With proper processing applied to DE images, calcified coronary arteries can be visualized, and data shows these images can be used in Agatston calcium scoring

-Calcium scoring done on proper-quality DE chest X-rays was considerably accurate, as confirmed by subsequent CT calcium scoring^[1]

-Using DE chest X-rays to engage in calcium scoring can expand coronary disease screening programs, keep dose low, and enhance patient care



[1] Yingnan Song, Hao Wu, Di Wen, Bo Zhu, Philipp Graner, Leslie Ciancibello, Haran Rajeswaran, Karma Salem, Mehrdad Hajmomenian, Robert C. Gilkeson and David L. Wilson, "Detection of coronary calcifications with dual energy chest X-rays: clinical evaluation," Int J Cardiovasc Imaging (2020).



Single Exposure Dual Energy Enables New Opportunities For X-ray Imaging



More confident, accurate image interpretation






Save time and benefit patients in Radiology, Emergency and Critical Care departments



No additional procedure or dose – same as chest X-ray



Leverage existing X-ray equipment and clinical techniques

 www.kaimaging.com
 (226) 215-9897
 info@kaimaging.com

Jay Potipcoe
Applications Specialist
jpotipcoe@kaimaging.com

Karim Karim
CTO
kkarim@kaimaging.com

ICEwebinars
A COOL SERIES FOR HOT TOPICS

NGUYEN VAN CONG

Has Attended the Following Webinar

**Improving Department Results with Single Exposure
Dual-Energy Subtraction X-Ray: Clinical, Operational
and Financial Benefits**

on

7/27/2021

CERTIFICATE OF ATTENDANCE

Sponsored by KA Imaging. Presented by Dr. Karim Karim, Ph.D., PEng, MBA and Jay Potipcoe, BA, MRT. Approved for 1 ARRT Category A CE credit by AHRA (Ref. LEC11337). Approved for 1 CRA credit. Expires 7/28/2022.



A handwritten signature in black ink, appearing to read 'Jayme Potipcoe'.



31°C Mưa nhỏ