
## GAAD score

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1. Henry Chan et al. Performance Evaluation of the Elecsys® GAAD Assay for the Detection of Hepatocellular Carcinoma Across Different Disease Stages and Etiologies - Presented at ISHVLD GHS 2021, 18-20 June, 2021. Global Hepatitis Summit.

\*navify Algorithm Suite is a Roche build software application hosted on navify platform which is designed to integrate, enable and consolidate individual (or group of) clinical high medical value algorithms, all in one unique access point on the LIS (Laboratory Information System), HIS (Hospital Information System) and/or EMR (Electronic Medical Record).

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**The clinical utility of Elecsys GAAD score in the diagnosis of hepatocellular carcinoma**


Chung-Feng Huang<sup>1</sup>, Ashish Sharma<sup>2</sup>, Ming-Lung Yu<sup>1</sup>  
<sup>1</sup>Hepatobiliary Division, Department of Internal Medicine, Kaohsiung Medical University Hospital, Kaohsiung Medical University, Kaohsiung, Taiwan  
<sup>2</sup>Clinical Development and Medical Affairs, Roche Diagnostics International AG, Rotkreuz, Switzerland

2022 APASL Single Topic Conference on Hepatocellular Carcinoma

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<b>Xơ gan, n (%)</b>	<b>32 (60.4)</b>
Child-Pugh A, n (%)	25 (47.2)
Child-Pugh B, n (%)	7 (13.2)
<b>Giai đoạn BCLC</b>	
BCLC 0, n (%)	7 (13.2)
BCLC A, n (%)	19 (35.8)
BCLC B, n (%)	12 (22.6)
BCLC C, n (%)	15 (28.3)
Kích thước khối u lớn nhất, cm (mean±SD)	5.0±3.6
AFP, ng/ml, (mean±SD)	2808±10788
<b>AFP &gt; 20 ng/ml, n (%)</b>	<b>23 (43.4)</b>
PIVKA-II, ng/mL, (mean±SD)	898±2077
<b>PIVKA &gt; 28.4 ng/mL, n (%)</b>	<b>37 (69.8)</b>
GAAD score, (mean±SD)	6.06±3.71
<b>GAAD score &gt; 2.57, n (%)</b>	<b>41 (77.4)</b>

## Tất cả ca

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
Điểm cắt : AFP-20 ng/mL, PIVKA-II 28.4 ng/mL, GAAD 2.57

PSAS

The clinical utility of Elecsys GAAD score in the diagnosis of hepatocellular carcinoma

2022 APASL Single Topic Conference on Hepatocellular Carcinoma

Chung-Feng Huang<sup>1</sup>, Ashish Sharma<sup>2</sup>, Ming-Lung Yu<sup>1</sup>  
<sup>1</sup>Hepatobiliary Division, Department of Internal Medicine, Kaohsiung Medical University Hospital, Kaohsiung Medical University, Kaohsiung, Taiwan  
<sup>2</sup>Clinical Development and Medical Affairs, Roche Diagnostics International AG, Rotkreuz, Switzerland

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
Huang CF et al. 2022 APASL STC on HCC, P161

## Thuật toán GAAD/ GALAD Training Study [Panel B]

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